

## Overview

---

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Interchangeable spindle spray arm with five direct-injection spindles for cleaning of small narrow-neck labware. For use in Labconco Upper or Lower Spindle Racks manufactured after October 2020.



**Catalog Number:** 4575800

## Specifications

---


- **Weight:** 2.0 lbs
- **Weight metric:** 0.9 kg
- **Dimensions:** 8.6" w x 0.5" d x 5.0" h
- **Dimensions metric:** 21.8 x 1.3 x 12.7 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower

## For Product Assistance

---



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative

