

## PRODUCT MODEL

Catalog # 8621110

# CentriVap Cold Trap -105 °C 230V 60 Hz US with Stainless Steel Lid



## All Features & Specifications

- Brushed, Stainless Steel and Powder-Coated Steel Exterior
- Illuminated Temperature Graph with Endpoint Indication (Operating temperature of -40 °C)
- HCFC-Free Refrigerant System
- 4-liter, Stainless Steel Collector
- Acrylic (-50 °C only) or Stainless Steel Lid with 2 Quick-Disconnect fittings for connection to a CentriVap Concentrator
- Stainless Steel Drain Port (with corrosion-resistant, nylon valve)
- Eco-cool natural refrigerant available on select models
- -50 °C models feature a 1/3 hp Compressor (-50 °C in less than 30 min.)
- -84 °C models feature two 1/3 hp Compressors (-84°C in less than 40 min.)
- -105 °C models feature two 1/3 hp Compressors (-105°C in less than 40 min.)

## Specifications

Weight	55.0
Dimensions	24.1" W x 23.9" D x 17.6" H
Metric Dimensions	61.2 W x 60.7 D x 44.7 H
Collector Temperature	-105° C
Lid Type	Stainless Steel
Plug Type	North America, 230V
Style	Benchtop
Type/application	Very Low Freezing Point Solvents
Conformance	UKCA, CE
Optional Accessories	Two-Port Manifold
Works With These Products	CentriVap Pro Vacuum Concentrators, CentriVap Standard Vacuum Concentrators
Eco-Cool	No