

SERO 12

SERO 12 VS SERO-FUGE BLOOD BANKING CENTRIFUGES


SERO 12

Sero Fuge™

Full height handles for easy to grip use*	Y	
Accepts Primary Tube Sets without inserts	Y	Y
Non-interfered by tube insert with slots	Y	
Available tubular mounts	Y	
Mounts	200	150
Height (in)	21	28
Weight (kg. lb)	3.26	3.65
Number of Tubes	4, 12	4, 12
Mounted Tube (inches)	4"	4"
Digital Display	Y	Y
Use Programmable Fixed Cycle	Y	Y
Max Spin Speed	3,000 RPM	3,000 RPM
Max Spin G-Force	1,500 g	1,500 g

* Tube at Standing

Uses All BD SeroFuge Rotors

