

# Test Tube Holder 14-20 mm



View online: <https://www.labconco.com/product/test-tube-holder-14-20-mm/8576>

## Overview



Test Tube Holders are add-on accessories to the Lower Test Tube Spindle Rack. The Holders keep tubes upright. Constructed of Type 304 stainless steel.

Test Tube Holder accommodates test tubes with 14-20 mm OD. Two may be stacked for increased capacity. Stacked holders can't accommodate an upper rack. Constructed of type 304 stainless steel. Requires Labconco Test Tube Spindle Conversion Kit 4575300 and Lower Spindle Rack manufactured after October 2020 (not included).

**Catalog Number:** 4574800


## Specifications

- **Weight:** 9.0 lbs
- **Weight metric:** 4.1 kg
- **Dimensions:** 8.4" w x 18.5" d x 5.3" h
- **Dimensions metric:** 21.3 x 47.0 x 13.5 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, MO 64132-2696  
Phone: (816) 333-8811

