

## Overview

---

Glassware Holders are made of Type 304 stainless steel or plastic and fit onto spindles to grip glassware and hold it in place.

Vinyl-coated. Type 304 stainless steel wire with height adjustment screw.



**Catalog Number:** 4424600P

## Specifications

---

- **Weight:** 0.4 lbs
- **Weight metric:** 0.2 kg
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower, Upper

## Description


---

## For Product Assistance

---



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, MO 64132-2696  
Phone: (816) 333-8811

