



Overview

Glassware Holders are made of Type 304 stainless steel or plastic and fit onto spindles to grip glassware and hold it in place.

Vinyl-coated. Type 304 stainless steel wire with height adjustment screw.

Catalog Number: 4424800P

Specifications


- **Weight:** 0.4 lbs
- **Weight metric:** 0.2 kg
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower, Upper

Description

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

