


Overview



Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Type 304 stainless steel. Capacity: 18 bottles. Grid size: 70 mm x 70 mm.
19.5" w x 8.8" d x 6.0" h (50 x 22 x 15 cm).

Catalog Number: 4589201


Specifications

- **Weight:** 2.0 lbs
- **Weight metric:** 0.9 kg
- **Dimensions:** 19.5"w x 8.8"d x 6"h
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower, Upper

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

