

Overview



Standard Open Racks and Spindle Racks for SteamScrubber, FlaskScrubber and FlaskScrubber Vantage Glassware Washers manufactured after October 2020.

Lower Small Spindle Rack has a high-capacity design with 60 spindles, 60 spindle plugs, 12 spindle arm plugs, spindle tool and wash arm for cleaning of small narrow-neck labware. Constructed of type 304 stainless steel. May be used in place of any lower rack in Labconco glassware washers manufactured after October 2020.

Catalog Number: 4668901


Specifications

- **Weight:** 21.0 lbs
- **Weight metric:** 9.5 kg
- **Dimensions:** 21.4" w x 20.9" d x 8.6" h
- **Dimensions metric:** 54.4 x 53.1 x 21.8 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Rack Location:** Lower

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

