

Overview



Standard Open Racks and Spindle Racks for SteamScrubber, FlaskScrubber and FlaskScrubber Vantage Glassware Washers manufactured after October 2020.

Lower Spindle Rack is designed for cleaning narrow-neck labware. Includes 36 spindles, 36 adjustable height clips, 18 glassware holders, 10 small spindles and clips, 36 spindle plugs, and 12 spindle arm plugs. Constructed of type 304 stainless steel. One is included with each FlaskScrubber and FlaskScrubber Vantage washer. For use in Labconco glassware washers manufactured after October 2020.

Catalog Number: 4668900


Specifications

- **Weight:** 21.0 lbs
- **Weight metric:** 9.5 kg
- **Dimensions:** 21.4" w x 20.9" d x 8.6" h
- **Dimensions metric:** 54.4 x 53.1 x 21.8 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Rack Location:** Lower

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

