

Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.



Mesh Basket is a full rack, open mesh basket insert that can be used for a variety of utensils and loose labware. Constructed of type 304 stainless steel. For use in Labconco Upper or Lower Standard Rack manufactured after October 2020.

Catalog Number: 4573400


Specifications

- **Weight:** 5.0 lbs
- **Weight metric:** 2.3 kg
- **Dimensions:** 18.6" w x 19.0" d x 6.7" h
- **Dimensions metric:** 47.2 x 48.3 x 17.0 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Standard
- **Rack Location:** Lower, Upper

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

