

Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

For Petri dishes, watch glasses and microtiter plates. Capacity: 24 plates. Grid size: 22 mm x 140 mm. Type 304 stainless steel. 19.2" x 8.6" x 5.9" h (49 x 22 x 15 cm).



Catalog Number: 4589701

Specifications

- **Weight:** 2.0 lbs
- **Weight metric:** 0.9 kg
- **Dimensions:** 19.2"w x 8.6"d x 5.9"h
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower, Upper

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

