

Overview



Standard Open Racks and Spindle Racks for SteamScrubber, FlaskScrubber and FlaskScrubber Vantage Glassware Washers manufactured after October 2020.

Upper Small Spindle Rack has a high-capacity design with 60 spindles, 60 spindle plugs, 12 spindle arm plugs, spindle tool and wash arm for cleaning of small narrow-neck labware. Constructed of type 304 stainless steel. For use in SteamScrubber washers manufactured after October 2020.

Catalog Number: 4668601


Specifications

- **Weight:** 21.0 lbs
- **Weight metric:** 9.5 kg
- **Dimensions:** 21.5" w x 20.8" d x 7.3" h
- **Dimensions metric:** 54.6 x 52.8 x 18.5 cm
- **Applicable Products:** SteamScrubber Glassware Washers
- **Rack Location:** Upper

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

