

# Upper Spindle Rack for FlaskScrubber and FlaskScrubber Vantage



View online: <https://www.labconco.com/product/upper-spindle-rack-for-flaskscrubber-and-flaskscrubber-vantage/8602>

## Overview



Standard Open Racks and Spindle Racks for SteamScrubber, FlaskScrubber and FlaskScrubber Vantage Glassware Washers manufactured after October 2020.

Upper Spindle Rack is designed for cleaning narrow-neck labware. Includes 36 spindles, 36 adjustable height clips, 18 glassware holders, 10 small spindles and clips, 36 spindle plugs, and 12 spindle arm plugs. Constructed of type 304 stainless steel. Includes upper slide rails. For use in FlaskScrubber and FlaskScrubber Vantage washers manufactured after October 2020.

**Catalog Number:** 4668500


## Specifications

- **Weight:** 23.0 lbs
- **Weight metric:** 10.4 kg
- **Dimensions:** 21.6" w x 20.8" d x 6.9" h
- **Dimensions metric:** 54.9 x 52.8 x 17.5 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers
- **Rack Location:** Upper

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, MO 64132-2696  
Phone: (816) 333-8811

