

# Upper Small Spindle Rack for FlaskScrubber and FlaskScrubber Vantage



View online: <https://www.labconco.com/product/upper-small-spindle-rack-for-flaskscrubber-and-flaskscrubber-vantage/8599>



## Overview

Standard Open Racks and Spindle Racks for SteamScrubber, FlaskScrubber and FlaskScrubber Vantage Glassware Washers manufactured after October 2020.

Upper Small Spindle Rack has a high-capacity design with 60 spindles, 60 spindle plugs, 12 spindle arm plugs, spindle tool and wash arm for cleaning of small narrow-neck labware. Constructed of type 304 stainless steel. Includes upper slide rails. For use in FlaskScrubber and FlaskScrubber Vantage washers manufactured after October 2020.

**Catalog Number:** 4668501


## Specifications

- **Weight:** 23.0 lbs
- **Weight metric:** 10.4 kg
- **Dimensions:** 21.5" w x 20.8" d x 7.3" h
- **Dimensions metric:** 54.6 x 52.8 x 18.5 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers
- **Rack Location:** Upper

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, MO 64132-2696  
Phone: (816) 333-8811

