

## Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Utensil Basket with cover securely holds miscellaneous labware such as stoppers, spatulas and other small lightweight items. Constructed of type 304 stainless steel. For use in Labconco Upper or Lower, Standard or Spindle Rack manufactured after October 2020.



**Catalog Number:** 4573100


## Specifications

- **Weight:** 3.0 lbs
- **Weight metric:** 1.4 kg
- **Dimensions:** 8.9" w x 8.9" d x 4.0" h
- **Dimensions metric:** 22.6 x 22.6 x 10.2 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower, Upper

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative

