

Overview



Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Utensil Holder for spatulas, spoons, scoops and other miscellaneous utensils. Utensils fit into six spacious compartments. Constructed of type 304 stainless steel. For use in Labconco Upper or Lower, Standard or Spindle Racks manufactured after October 2020.

Catalog Number: 4575200


Specifications

- **Weight:** 5.0 lbs
- **Weight metric:** 2.3 kg
- **Dimensions:** 8.2" w x 8.2" d x 8.0" h
- **Dimensions metric:** 20.8 x 20.8 x 20.3 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower, Upper

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

