

Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Vial Insert securely holds up to 48 each 15 x 45 mm vials during washing. Includes retainer lid with handle, support holder and six spindle arms. May use up to two per Spindle Rack at one time. Constructed of type 304 stainless steel. For use in Labconco Lower Spindle Racks manufactured after October 2020.



Catalog Number: 4574500


Specifications

- **Weight:** 8.0 lbs
- **Weight metric:** 3.6 kg
- **Dimensions:** 8.7" w x 19.2" d x 7.9" h
- **Dimensions metric:** 22.1 x 48.8 x 20.1 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

