

## Overview



Add mobility to freestanding SteamScrubber and FlaskScrubber Glassware Washers. Allow washers to be plugged into a 20 amp electrical wall outlet and connected to a standard laboratory sink faucet. Include four 2-inch diameter rubber casters, two aerator adapters for connection to standard 55/64"-27 male and 15/16"-27 female laboratory faucets, drain hose and installation instructions. Not for use with undercounter washers or FlaskScrubber Vantage Series Glassware Washers.

Freestanding-to-Mobile Conversion Kit adds mobility to freestanding washer models. Allows washers to be plugged into a 20A electrical wall outlet & connected to a standard laboratory sink faucet. For use with Labconco FlaskScrubber and SteamScrubber Washers manufactured after October 2020. Not for use with FlaskScrubber Vantage models. For 230V operation.

**Catalog Number:** 4596111


## Specifications

- **Weight:** 60.0 lbs
- **Weight metric:** 27.2 kg
- **Dimensions:** 15.0" w x 4.0" d x 4.0" h
- **Dimensions metric:** 38.1 x 10.2 x 10.2 cm
- **Electrical:** 230V
- **Applicable Products:** FlaskScrubber Glassware Washers, SteamScrubber Glassware Washers

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative

