

Overview

Heater Cover provides barrier layer between labware and heating coils in Labconco FlaskScrubber and FlaskScrubber Vantage washers manufactured after October 2020. Unique design provides additional filtration layer to protect from broken glass. Constructed of type 304 stainless steel.



Catalog Number: 4679301


Specifications

- **Weight:** 6.0 lbs
- **Weight metric:** 2.7 kg
- **Dimensions:** 21.4" w x 20.4" d x 0.75" h
- **Dimensions metric:** 54.4 x 51.8 x 1.9 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

