

## Overview

---



Heater Cover provides barrier layer between labware and heating coils in Labconco SteamScrubber washers manufactured after October 2020. Unique design provides additional filtration layer to protect from broken glass. Constructed of type 304 stainless steel.

**Catalog Number:** 4679300

## Specifications

---


- **Weight:** 6.0 lbs
- **Weight metric:** 2.7 kg
- **Dimensions:** 21.4" w x 20.4" d x 0.75" h
- **Dimensions metric:** 54.4 x 51.8 x 1.9 cm
- **Applicable Products:** SteamScrubber Glassware Washers

## For Product Assistance

---



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative

