



## Overview

Automatic dispenser kits and storage for detergent and neutralizing acid rinse with Labconco glassware washers.

Liquid Detergent and Rinse Aid Dispenser Kit connects to any SteamScrubber or FlaskScrubber washer to automate liquid detergent and neutralizing acid rinse aid dispensing. Includes peristaltic pump, float switch, tubing, one-gallon plastic container for detergent, and one-liter plastic container for rinse aid. Detergent and rinse aid not included. For use with washers manufactured after October 2020.

**Catalog Number:** 4679700


## Specifications

- **Weight:** 12.0 lbs
- **Weight metric:** 5.4 kg
- **Dimensions:** 9.25" w x 9.25" d x 11.8" h
- **Dimensions metric:** 23.5 x 23.5 x 30.0 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, SteamScrubber Glassware Washers

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative

