

## PRODUCT MODEL

Catalog # 780401000

# Lyph-Seal Tray Dryer 115V US



## All Features & Specifications

- Stoppering/ Sealing of a non-rigid container
- Ability to seal objects that were previously excluded from glass vial stoppering.
- Dry and Seal in the final container, without handling outside of a controlled environment.
- Compatibility with all Labconco Freeze Dryers equipped with LyoWorks OS (2017 to current)
- Includes 4 Heated shelves and sample temperature probe.
- Heated product shelves allow temperature set points between -20°C and +60°C in 1-degree increments.
- Heated product shelves allow for expanded control and optimization of the Lyophilization process.
- Shelves are independently controlled to a temperature set point, preventing shelf-to-shelf variability when thermal masses differ.
- 3-way valve allows for simple selection of “Drying, sealing, and venting” conditions within chamber.

## Specifications

---

<b>Dimensions</b>	15.6" W x 19.2" D x 17.9" H
<b>Electrical</b>	115 volts, 50/60 Hz, 4 amps
<b>Height</b>	0
<b>Shipping Length</b>	0
<b>Metric Dimensions</b>	39.6 W x 48.8 D x 45.5 H
<b>Metric Weight (kg)</b>	62.6
<b>Width</b>	0
<b>Plug Type</b>	North America, 115V
<b>Vacuum Pump Displacement Required</b>	98 L/min or Larger Displacement
<b>Applicable Products</b>	Freezone 12 Liter, Freezone 18 Liter, Freezone 2.5 Liter, Freezone 4.5 Liter, Freezone 6 Liter, Freezone 8 Liter
<b>Conformance</b>	CAN/CSA C22.2

---