

MacroTube® VIVID™

- **Vivid Printed Graduations**
The only printed 5mL tube!
- **Resealable Sterile Packaging**
- **RNase, DNase, Pyrogen-free**
- **Ultra-Clear Polypropylene**
- **Snap-Apart Racks**
2 x 25 tubes - 50 total

Specifications

Material (Tube):	Polypropylene
Material (Cap):	High Density Polyethylene
Maximum RFC:	25,000xg
Temp Range:	-86°C to 100°C
Dimensions: (Without cap & rim)	16 x 54mm
Dimensions: (Withcap & rim)	20 x 64mm

DNA, RNase, DNase, & Pyrogen Free



**Printed
Graduations**

The only printed 5mL tube!

**5mL
Centrifuge
Tubes**

**Ultra
Clear**

Snap Apart Racks
2x25 Tubes per Rack



Cat. No.	MacroTube® Vivid - 5mL Tubes w/ Vivid Printed Graduations	QTY
C2545	5mL Macro Centrifuge Tube, White Printed Graduations and Marking Area, with Screw Caps Attached, Sterile, 10 Snap-Apart Racks of 50 Tubes	500/cs

