

# Powerful at every turn

## Save up to 30% with Multifuge Pro and Megafuge Plus General Purpose Benchtop Centrifuge Packages

### Thermo Scientific™ Multifuge™ Pro and Megafuge™ Plus General Purpose Centrifuges

- Available in 1.6 L benchtop, 4 L benchtop, or floor standing options
- Touchscreen—easy to use, helps save time, and promotes sample security with secure access and logging capabilities
- Improved ergonomics—loading height decreased and less force needed to close the lid
- Thermo Scientific™ AUTO-LOCK™ II Rotor Exchange system—switch applications and rotors safely with no tools required
- Spin up to 40 x 50 mL conical tubes, 196 x 5 to 7 mL blood tubes and 24 microplates in Thermo Scientific™ TX-1000 rotor (4 L units) and 16 x 50 mL conical tubes and 84 x 5 to 7 mL and 10 microplates in Thermo Scientific™ H-Flex 1 rotor (1.6 L units)
- Certified biocontainment—Thermo Scientific™ ClickSeal™ Biocontainment Lids provide one-handed sample biocontainment



The family of Multifuge Pro and Megafuge Plus General Purpose Centrifuges includes the 1.6 L & 4 L benchtop centrifuges and floor standing options.




Auto-Lock Rotor Exchange



ClickSeal Biocontainment Lids



Fiberlite Carbon Fiber Rotors



Unmatched capacity and application flexibility

# Multifuge X1 Pro and Megafuge ST1 Plus Centrifuges



**Multifuge X1R Pro Centrifuge**



**Megafuge ST1 Plus Centrifuge**

## Swing-out rotor packages

Thermo Scientific™ TX-400 blood tube and cell culture packages include the following with one part number for simple ordering

- Your choice of centrifuge
- TX-400 or H-Flex 1 rotor and buckets with ClickSeal biocontainment lids
- Adapters for application

## Maximized capacity for your applications:

- Cell culture and molecular biology applications: spin up to 16 x 50 mL conical tubes or 10 standard microplates
- Clinical applications: spin up to 84 x 5/7 mL blood tubes

Centrifuge	Cat. No.
<b>1.6 L high-performance cell culture packages</b>	
<b>TX-400 Rotor: 16 x 50 mL conical tubes up to 5,000 rpm/4,696 x g</b>	
Multifuge X1R Pro Refrigerated TX-400 Cell Culture Package	75016008
Multifuge X1 Pro Ventilated TX-400 Cell Culture Package	75016024US
Megafuge ST1R Plus Refrigerated TX-400 Cell Culture Package	75016010
Megafuge ST1 Plus Ventilated TX-400 Cell Culture Package	75016026US
<b>1.6 L high-performance blood tube packages</b>	
<b>TX-400 Rotor: 76 x 5/7 mL blood tubes up to 5,000 rpm/4,696 x g</b>	
Megafuge ST1R Plus Refrigerated TX-400 Blood Tube Package	75016002
Megafuge ST1 Plus Ventilated TX-400 Blood Tube Package	75016003
<b>H-Flex 1 cell culture - 16 x 50 mL conical tubes and 10 microplates up to 4,700 rpm/4,297 x g</b>	
Multifuge X1R Pro Refrigerated H-Flex Cell Culture Package	75017008
Multifuge X1 Pro Ventilated H-Flex Cell Culture Package	75017024US
<b>H-Flex 1 blood tube - 84 x 5/7 mL blood tubes and 10 microplates up to 4,700 rpm/4,297 x g</b>	
Megafuge ST1R Plus Refrigerated H-Flex Blood Tube Package	75017002
Megafuge ST1 Plus Ventilated H-Flex Blood Tube Package	75017003

### Package includes



- TX-400 cell culture rotor package includes:
- Your choice of centrifuge
  - TX-400 rotor cross
  - TX-400 buckets (set of 4)
  - TX-400 ClickSeal Biocontainment Lids (set of 4)
  - Adapters 16 x 50 mL conical (set of 4)
  - Adapters 36 x 15 mL conical (set of 4)



- TX-400 blood tube rotor package includes:
- Your choice of centrifuge
  - TX-400 rotor cross
  - TX-400 buckets (set of 4)
  - TX-400 ClickSeal Biocontainment lids (set of 4)
  - Adapters 56 x 10 mL (set of 4)
  - Adapters 76 x 5/7 mL (set of 4)



- H-FLEX 1 package includes:
- Your choice of centrifuge
  - H-Flex 1 rotor cross, set of 2 buckets and set of 2 ClickSeal Biocontainment lids
- Cell culture and microplate package includes:
  - Adapters 16 x 50 mL conical (set of 2)
  - Adapters 36 x 15 mL conical (set of 2)
  - Microplate carriers (set of 2)
- Blood processing and microplate package includes:
  - Adapters 66 x 10 mL (set of 2)
  - Adapters 84 x 5/7 mL (set of 2)
  - Microplate carriers (set of 2)

# Multifuge X4 Pro and Megafuge ST4 Plus benchtop and floor-model centrifuges

## Swing-out rotor packages

Our TX-1000 blood tube and cell culture packages include the following with simple, one part number ordering

- Your choice of 4 L benchtop or floor standing centrifuge
- TX-1000 Rotor
- TX-1000 Rotor Buckets
- ClickSeal Biocontainment Lids
- Adapters for application

Our TX-1000 microplate and TX-750 blood tube and cell culture packages are available, but the centrifuge must be purchased separately.



Megafuge ST4R Plus Centrifuge



Multifuge X4R Pro Centrifuge

## Maximized capacity for your applications:

- Cell culture and molecular biology applications: spin up to 40 x 50 mL conical tubes or 24 standard microplates
- Clinical applications: spin up to 196 x 5 to 7 mL blood tubes

Centrifuge	Cat. No.	Package includes	
<b>4 L high-performance and high-throughput cell culture packages</b>			
<b>40 x 50 mL tubes up to 4,200 rpm/4,122 x g (refrigerated) or up to 3,800 rpm/3,374 x g (ventilated)</b>			
Multifuge X4R Pro Refrigerated TX-1000	75016045		
Multifuge X4 Pro Ventilated TX-1000	75016053		
Megafuge ST4R Plus Refrigerated TX-1000	75016049		
Megafuge ST4 Plus Ventilated TX-1000	75016050		
Multifuge X4FR Pro Floor Standing Refrigerated TX-1000	75016047		
Multifuge X4F Pro Floor Standing Ventilated TX-1000	75016048		
Megafuge ST4FR Plus Floor Standing Refrigerated TX-1000	75017946		
Megafuge ST4 Plus Floor Standing Ventilated TX-1000	75017960		
<b>4 L high performance and high throughput blood tube packages</b>			
<b>196 x 5/7 ml or 148 x 10 mL blood tubes up to 4,200 rpm/4,122 x g (refrigerated) or up to 3,800 rpm/3374 x g (ventilated)</b>			
Multifuge X4R Pro Refrigerated Blood Tube Package	75016033		
Multifuge X4 Pro Ventilated TX-1000 Blood Tube Package	75016034		
Megafuge ST4R Plus Refrigerated TX-1000 Blood Tube Package	75016037		
Megafuge ST4 Plus Ventilated TX-1000 Blood Tube Package	75016038		
Megafuge ST4FR Plus Floor Standing Refrigerated TX-1000 Blood Tube Package	75016059		
Megafuge ST4F Plus Floor Standing Ventilated TX-1000 Blood Tube Package	75017964		
<b>Add-on rotor packages—add flexibility to your centrifuge</b>			
<b>Purchase your centrifuge and add on a package below</b>			
TX-750 Cell Culture Package	75016110		
TX-750 Blood Tube Package	75016109		

TX-1000 Microplate Package	75016104
----------------------------	----------

- TX-1000 cell culture rotor packages include:
- Your choice of centrifuge
  - TX-1000 rotor cross
  - TX-1000 buckets (set of 4)
  - TX-1000 ClickSeal Biocontainment lids (set of 4)
  - Adapters 40 x 50 mL conical (set of 4)
  - Adapters 96 x 15 mL conical (set of 4)

- TX-1000 blood tube rotor packages includes:
- Your choice of centrifuge
  - TX-1000 rotor cross
  - TX-1000 buckets (set of 4)
  - TX-1000 ClickSeal Biocontainment Lids (set of 4)
  - Adapters 148 x 10 mL (set of 4)
  - Adapters 196 x 5 to 7 mL (set of 4)

- TX-750 packages includes:
- TX-750 rotor cross
  - TX-750 buckets (set of 4)
  - TX-750 ClickSeal Biocontainment Lids (set of 4)
  - TX-750 50 mL adapter (set of 4)
  - TX-750 15 mL adapter (set of 4)
- OR
- TX-750 10 mL blood tube (set of 4)
- TX-750 5 to 7 mL blood tube (set of 4)

- TX-1000 package includes:
- TX-1000 rotor cross
  - TX-1000 buckets (set of 4)
  - TX-1000 ClickSeal Biocontainment Lids (set of 4)
  - Microplate adapters (set of 4)

## Increase your package flexibility even further with add-ons

After you choose your application-focused package, simply select from complementary Thermo Scientific™ products to further enhance your lab's processing capabilities.



### Add-on rotors and adapters

Thermo Scientific™ rotors and adapters	Cat. No.
Fiberlite F14-6 x 250 LE rotor	75773662
Fiberlite F13-14 x 50cy rotor	75773661
Fiberlite F15-6 x 100y rotor	75773698
Fiberlite F15-8 x 50cy rotor	75773663
HIGHConic II rotor	75773620
Microliter 48 x 2 rotor	75773602
Microliter 30 x 2 rotor	75773652
MicroClick 30 x 2 rotor (prep kit compatible)	75775719
TX-1000 microplate adapters*	75777303
M-20 rotor with standard microplate carriers	75773624
Sealed microplate carriers for M-20 rotor*	75773625
MicroClick 18 x 5 rotor	75775765

\*Rotor sold separately. Check on rotor and accessory compatibility with selected centrifuge.



Thermo Scientific™ mySPIN™ 12 and 6 mini centrifuges

### Add-on centrifuges

Centrifuges	Includes	Max speed	Max RCF	Cat. No.
mySPIN 12 mini centrifuge	<ul style="list-style-type: none"> <li>• 12 x 1.5/2.0 mL tube rotor</li> <li>• 32 x 0.2 mL tube rotor (singles or 4 x 8 strip)</li> <li>• 0.2 mL and 0.5 mL tube adapters (sets of 12)</li> </ul>	12,500 rpm	9,800 x g	75774081
mySPIN 6 mini centrifuge	<ul style="list-style-type: none"> <li>• 6 x 1.5/2.0 mL tube rotor</li> <li>• 16 x 0.2 mL tube rotor (singles or 2 x 8 strip)</li> <li>• 0.2 mL and 0.5 mL tube adapters (sets of 6)</li> </ul>	6,000 rpm	2,000 x g	75774061

Find out more at [thermofisher.com/acceleratingscience](https://www.thermofisher.com/acceleratingscience)