

Overview

Automatic dispenser kits and storage for detergent and neutralizing acid rinse with Labconco glassware washers.

Acid Rinse Dispenser Kit connects to any SteamScrubber or FlaskScrubber washer to automate neutralizing acid rinse aid dispensing. Includes peristaltic pump, float switch, tubing, and one-liter plastic container. Rinse aid not included. For use with washers manufactured after October 2020.



Catalog Number: 4595120


Specifications

- **Weight:** 3.0 lbs
- **Weight metric:** 1.4 kg
- **Dimensions:** 3.25" w x 3.25" d x 8.0" h
- **Dimensions metric:** 8.3 x 8.3 x 20.3 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, SteamScrubber Glassware Washers

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

