

Picus[®] 2 Electronic Pipette

The Next Generation of Pipetting

Picus[®] 2 Electronic Pipette


- Reliable performance and reproducible results
- Easily accessible pipetting modes and customizable programs speed up repetitive work
- Lightweight and ergonomic design reduces muscle strain and provides comfortable pipetting
- Sartorius pipetting mobile app provides configurable protocols and pipette management
- Open connectivity feature provides integration of Picus[®] 2 pipette into existing lab operating systems
- Service and calibration reminders
- Password protection and distinct administration user profiles protect against unauthorized modifications
- Calibration adjustment to provide accurate results even with the most challenging liquid
- Easy to clean and decontaminate; lower parts of the pipette can be decontaminated by autoclaving (excluding 1,200 µL multichannel models), with UV, or with disinfectants
- Battery typically lasts a full working day and can be fully recharged in an hour

Sartorius Mobile App:

- Sample preparation workflows and pipette management
- Step-by-step workflow guidance and automatic selection of correct settings for Picus[®] 2
- Pipette updates
- Ever-growing database of configurable protocols



Learn More



Ordering Information: Picus® 2


Order Codes	Channels	Volume Range (µl)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace Filter Tips (µl)
LH-747021	1	0.2-10	-	-	10	10
LH-747041	1	5-120	721008	721018	200	120
LH-747061	1	10-300	721007	721017	350	300
LH-747081	1	50-1,000	721006	721016	1,000	1,000
LH-747101	1	100-5,000	721005	721015	5,000	5,000
LH-747111	1	500-10,000	721005	721015	10,000	-
LH-747321	8	0.2-10	-	-	10	10
LH-747341	8	5-120	721008	721018	200	120
LH-747361	8	10-300	721007	721017	350	300
LH-747391	8	50-1,200	721006	721016	1,200	1,200
LH-747421	12	0.2-10	-	-	10	10
LH-747441	12	5-120	721008	721018	200	120
LH-747461	12	10-300	721007	721017	350	300
LH-747491	12	50-1,200	721006	721016	1,200	1,200



Pipette Stands and Charging Options

Order Code	Name
LH-735001	USB charging adapter
730981	Charging stand for 1 Picus® pipette
730991	Charging carousel for 4 Picus® pipettes
725620	Linear Stand (non-charging)

Download the App



For more information, please contact your Sartorius Sales Representative: