

special offer
SUMMER
SALE

SHOP NOW
 up to 50 % off



qTOWER³ Real-Time PCR Systems

High Performance DNA Amplification and Quantification

- Patented fiber optic system for ideal real-time polymerase chain reaction signals, market leading temperature control, capable of up to six-fold channel multiplexing
- Open platform accommodates numerous plasticware and kits/mastermixes, **10-year** warranty on optical components; **3-year** warranty on the remainder of the system and completely maintenance and calibration free
- Solid silver sample block provides heating/cooling rates of up to 8.0°C/sec and 6.0°C/sec (respectively), and temperature homogeneity of ±0.15°C (96-well models only)
- License-free control and analysis software allows for installation on multiple computers



Ordering Information

| Order Number | Description | List Price | Promo Price |
|--------------|---|-------------|-------------|
| SP-1037 | qTOWER3 G: 96-well qPCR system, incl. 2 channels and laptop computer | \$32,649.00 | \$17,995.00 |
| SP-1038 | qTOWER3 G: 96-well qPCR system, incl. 6 channels and laptop computer | \$51,960.00 | \$23,995.00 |
| SP-1039 | qTOWER384 G: 384-well qPCR system, incl. 2 channels and laptop computer | \$54,556.00 | \$28,995.00 |
| SP-1040 | qTOWER384 G: 384-well qPCR system, incl. 6 channels and laptop computer | \$73,867.00 | \$34,995.00 |

Promotions valid only in North America. Orders must be placed by 9/30/2021.

For full product specifications, brochures and application notes, please visit our website: www.analytik-jena.us