



## Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Full-coverage spray nozzle distributes water to maximize cleanability of wide-mouth labware. Replaces standard spindle in Labconco Upper or Lower Spindle Racks manufactured after October 2020 (not included). Includes one spray nozzle and one spindle.

**Catalog Number:** 4575700


## Specifications

- **Weight:** 2.0 lbs
- **Weight metric:** 0.9 kg
- **Dimensions:** 0.13" w x 0.13" d x 3.7" h
- **Dimensions metric:** 0.3 x 0.3 x 9.4 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower, Upper

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative

