

Model

MLR102 Medical-grade Countertop Refrigerator


Specifications

Application, Rating and Electrical Data		
Application	Medication, vaccine, and other medical and scientific products	
Model Type	Countertop	
Storage Volume	1.8 cu ft (51 L)	
Temperature Range	2° to 8°C	
Electrical Power	115V 60Hz	
Maximum Current	0.75 Amps	
Building Supply Rating	15 amp dedicated ground circuit	
Power Plug / Power Cord Length	NEMA 5-15 hospital grade plug 6 ft (1.8 m)	
Certification / Agency Listing	UL/cUL Certified	
Indoor/Outdoor Use	Indoor use only	
Application Environment	Non-corrosive, non-flammable, non-explosive, 10°C to 32°C (50°F to 90°F)	
Refrigeration System		
Cooling System	AdvanceCore [™] solid-state thermoelectric heat pump technology with forced-air circulation	
Refrigerant Charge	R744 (CO ₂), non-flammable, zero ozone depleting, low toxicity	
Performance Data		
Uniformity	+/-1.0°C	
Noise Level	40dB	
Energy Requirements	0.62 kWh/24hrs @ 4°C	
Warranty		
Warranty (US and Canada)	2 years parts and labor	


	Controller	
	Interface/Display	Color LCD
	Controller Type	Microprocessor with menu-based input and display. Includes USB port for data retrieval.
	Security	Lockable door, password protected calibration settings
	Communication Ports	USB Port
	High / Low Alarms	Fully adjustable
	High / Low Alarm Test	Manual
	Door Ajar Alarm	Yes, 1 min
	Power Failure Alert	Yes
	Event and Temperature Log	All alarms, download via USB
	Min/Max Temperature Display	Yes, in °C/°F
	Temperature Units	°C/°F
	Battery Back-up	Battery back-up for monitoring and data logging (CR123A)
Dimensions and Construction		1
	Interior (w x h x d)	14.3 x 18.4 x 12.2 in 363 x 467 x 309mm
	Exterior (w x h x d) (Includes leveling feet)	20.5 x 23 x 22.1 in 520 x 584 x 561 mm
	Insulation	Foamed
	Exterior / Interior Finish	Bacteria-resistant powder coating - exterior Aluminum - interior
	Access Ports	Rear access port
	Door	Solid, left-hinged
	Interior Storage /	3 stainless steel shelves (2 adjustable) 13.5" x 10.2" (343 x 259mm)
	Capacity	40 lbs (18 kg) max capacity / shelf
	Leveling Feet	
		40 lbs (18 kg) max capacity / shelf
	Leveling Feet	40 lbs (18 kg) max capacity / shelf Yes
	Leveling Feet Net Weight	40 lbs (18 kg) max capacity / shelf Yes 75 lb (34 kg)









Note: Height with feet retracted to 0.75" (19mm)







Helmer Scientific | 14400 Bergen Boulevard | Noblesville, IN 46060 | USA +1.317.773.9073 | www.helmerinc.com

