



Centrifuges

Centrifugation packages

Compact and general purpose centrifuges

Discover our Thermo Scientific™ Centrifuge Packages designed to support cell culture, molecular biology, blood/clinical and multi-discipline applications, helping to save you time by combining our general purpose centrifuges

with the necessary components to run specific protocols in one convenient bundle. Select different capacities in our compact centrifuge and general purpose centrifuges with different rotor and adapter options.

Cell culture applications



Thermo Scientific™ Megafuge™ and Multifuge™ Packages



Image: Megafuge 8 centrifuge

Maximum capacities

4x
50 mL

8x
15 mL

Package	Centrifuge	Accessories included	Cat. No.
Megafuge 8 Ventilated Cell Culture Package	Megafuge 8 ventilated centrifuge, 230 V	- TX-150 rotor (75005701) - TX-150 round buckets (set of 4) (75005702) - TX-150 ClickSeal biocontainment lids (set of 4) (75005707)	50140758
Megafuge 8R Refrigerated Cell Culture Package	Megafuge 8R refrigerated centrifuge, 230 V	- Adapters 50 mL conical (set of 4) (75005736) - Adapters 2 x 15 mL conical (set of 4) (75005737)	50144363



Image: Megafuge 8 centrifuge

8x
50 mL

8x
15 mL

Package	Centrifuge	Accessories included	Cat. No.
Megafuge 8 Ventilated High-Capacity Cell Culture Package	Megafuge 8 ventilated centrifuge, 230 V	- TX-150 rotor (75005701) - TX-150 conical buckets (set of 4) (75005703) - Adapters 15 mL (set of 8) (75005808)	50147012



Image: Megafuge ST1 Plus centrifuge

16x
50 mL


36x
15 mL

Package	Centrifuge	Accessories included	Cat. No.
Multifuge X1 Pro Ventilated TX-400 Cell Culture Package	Multifuge X1 Pro ventilated centrifuge, 230 V	- TX-400 rotor cross - TX-400 buckets (set of 4) - TX-400 ClickSeal biocontainment lids (set of 4)	75016024
Multifuge X1R Pro Refrigerated TX-400 Cell Culture Package	Multifuge X1R Pro refrigerated centrifuge, 230 V		75016023
Megafuge ST1 Plus Ventilated TX-400 Cell Culture Package	Megafuge ST1 Plus ventilated centrifuge, 230 V	- Adapters for 4 x 50 mL conical tubes (set of 4) - Adapters for 9 x 15 mL conical tubes (set of 4)	75016026
Megafuge ST1R Plus Refrigerated TX-400 Cell Culture Package	Megafuge ST1R Plus refrigerated centrifuge, 230 V		75016025

Thermo Scientific™ Centrifuge Packages



Image: X4T centrifuge



28x
50 mL

56x
15 mL

Package	Centrifuge	Accessories included	Cat. No.
X4T TX-750 Cell Culture Package (GUO)	X4T ventilated centrifuge, GUO (general use only), 230 V	- TX-750 rotor cross - TX-750 bucket (set of 4) - TX-750 ClickSeal biocontainment lids (set of 4) - Adapters 7 x 50 mL conical (set of 4) - Adapters 14 x 15 mL conical (set of 4)	75015816
X4T TX-750 Cell Culture Package (IVD-MD)	X4T ventilated centrifuge, IVD-MD, 230 V		75015812
X4TR TX-750 Cell Culture Package (GUO)	X4TR refrigerated centrifuge, GUO, 230 V		75015848
X4TR TX-750 Cell Culture Package (IVD-MD)	X4TR refrigerated centrifuge, IVD-MD, 230 V		75015844
S4T TX-750 Cell Culture Package (GUO)	S4T ventilated centrifuge, GUO, 230 V		75015806
S4T TX-750 Cell Culture Package (IVD-MD)	S4T ventilated centrifuge, IVD-MD, 230 V		75015802
S4TR TX-750 Cell Culture Package (GUO)	S4TR refrigerated centrifuge, GUO, 230 V		75015832
S4TR TX-750 Cell Culture Package (IVD-MD)	S4TR refrigerated centrifuge, IVD-MD, 230 V		75015828



Image: X4T centrifuge



40x
50 mL

96x
15 mL

	Package	Centrifuge	Accessories included	Cat. No.
Benchtop	X4T TX-1000 Cell Culture Package (GUO)	X4T ventilated centrifuge, GUO, 230 V	- TX-1000 rotor cross - TX-1000 buckets (set of 4) - TX-1000 ClickSeal biocontainment lids (set of 4) - Adapters 10 x 50 mL conical (set of 4) - Adapters 24 x 15 mL conical (set of 4)	75015814
	X4T TX-1000 Cell Culture Package (IVD-MD)	X4T ventilated centrifuge, IVD-MD, 230 V		75015810
	X4TR TX-1000 Cell Culture Package (GUO)	X4TR refrigerated centrifuge, GUO, 230 V		75015846
	X4TR TX-1000 Cell Culture Package (IVD-MD)	X4TR refrigerated centrifuge, IVD-MD, 230 V		75015842
	S4T TX-1000 Cell Culture Package (GUO)	S4T ventilated centrifuge, GUO, 230 V		75015804
	S4T TX-1000 Cell Culture Package (IVD-MD)	S4T ventilated centrifuge, IVD-MD, 230 V		75015800
	S4TR TX-1000 Cell Culture Package (GUO)	S4TR refrigerated centrifuge, GUO, 230 V		75015830
	S4TR TX-1000 Cell Culture Package (IVD-MD)	S4TR refrigerated centrifuge, IVD-MD, 230 V		75015826
Floor	X4FR TX-1000 Cell Culture Package (GUO)	X4FR refrigerated floor-standing centrifuge, GUO, 230 V		75015838
	X4FR TX-1000 Cell Culture Package (IVD-MD)	X4FR refrigerated floor-standing centrifuge, IVD-MD, 230 V		75015836
	S4FR TX-1000 Cell Culture Package (GUO)	S4FR refrigerated floor-standing centrifuge, GUO, 230 V		75015822
	S4FR TX-1000 Cell Culture Package (IVD-MD)	S4FR refrigerated floor-standing centrifuge, IVD-MD, 230 V		75015820

Clinical applications



Thermo Scientific™ Megafuge™ Centrifuge Packages



Image: Megafuge 8 centrifuge

Maximum capacities

24x
5/7 mL

16x
10 mL

Package	Centrifuge	Accessories included	Cat. No.
Megafuge 8 Ventilated Routine Clinical Package	Megafuge 8 ventilated centrifuge, 230 V	- TX-150 rotor (75005701) - TX-150 round buckets (set of 4) (75005702) - TX-150 ClickSeal biocontainment lids (set of 4) (75005707)	50140775
Megafuge 8R Refrigerated Routine Clinical Package	Megafuge 8R refrigerated centrifuge, 230 V	- Adapters 5/7 mL blood collection (set of 4) (75005739) - Adapters 10/12 mL blood collection or 13 mL urine (set of 4) (75003504)	50144445



Image: Megafuge ST1 Plus centrifuge

76x
5/7 mL

56x
10 mL

Package	Centrifuge	Accessories included	Cat. No.
Megafuge ST1 Plus Ventilated TX-400 Blood Tube Package	Megafuge ST1 Plus-MD ventilated centrifuge, 230 V	- TX-400 rotor cross - TX-400 buckets (set of 4) - TX-400 biocontainment lids (set of 4)	75016003
Megafuge ST1R Plus Refrigerated TX-400 Blood Tube Package	Megafuge ST1R Plus-MD refrigerated centrifuge, 230 V	- Adapters 19 x 5/7 mL blood tubes (set of 4) - Adapters 14 x 10 mL blood tubes (set of 4)	75016019

Thermo Scientific™ Centrifuge Packages



Image: S4T centrifuge




108x
5/7 mL

84x
10 mL

Package	Centrifuge	Accessories included	Cat. No.
X4T TX-750 Blood Tube Package (GUO)	X4T ventilated centrifuge, GUO, 230 V	- TX-750 rotor cross - TX-750 buckets (set of 4) - TX-750 ClickSeal biocontainment lids (set of 4) - Adapters 21 x 10 mL blood tubes (set of 4) - Adapters 27 x 5/7 mL blood tubes (set of 4)	75015817
X4T TX-750 Blood Tube Package (IVD-MD)	X4T ventilated centrifuge, IVD-MD, 230 V		75015813
X4TR TX-750 Blood Tube Package (GUO)	X4TR refrigerated centrifuge, GUO, 230 V		75015849
X4TR TX-750 Blood Tube Package (IVD-MD)	X4TR refrigerated centrifuge, IVD-MD, 230 V		75015845
S4T TX-750 Blood Tube Package (GUO)	S4T ventilated centrifuge, GUO, 230 V		75015807
S4T TX-750 Blood Tube Package (IVD-MD)	S4T ventilated centrifuge, IVD-MD, 230 V		75015803
S4TR TX-750 Blood Tube Package (GUO)	S4TR refrigerated centrifuge, GUO, 230 V		75015833
S4TR TX-750 Blood Tube Package (IVD-MD)	S4TR refrigerated centrifuge, IVD-MD, 230 V		75015829



Image: S4T centrifuge



196x
5/7 mL

148x
10 mL

	Package	Centrifuge	Accessories included	Cat. No.
Benchtop	X4T TX-1000 Blood Tube Package (GUO)	X4T ventilated centrifuge, GUO, 230 V	- TX-1000 rotor cross - TX-1000 buckets (set of 4) - TX-1000 ClickSeal biocontainment lids (set of 4) - Adapters 49 x 5/7 mL blood tubes (set of 4) - Adapters 37 x 10 mL blood tubes (set of 4)	75015815
	X4T TX-1000 Blood Tube Package (IVD-MD)	X4T ventilated centrifuge, IVD-MD, 230 V		75015811
	X4TR TX-1000 Blood Tube Package (GUO)	X4TR refrigerated centrifuge, GUO, 230 V		75015847
	X4TR TX-1000 Blood Tube Package (IVD-MD)	X4TR refrigerated centrifuge, IVD-MD, 230 V		75015843
	S4T TX-1000 Blood Tube Package (GUO)	S4T ventilated centrifuge, GUO, 230 V		75015805
	S4T TX-1000 Blood Tube Package (IVD-MD)	S4T ventilated centrifuge, IVD-MD, 230 V		75015801
	S4TR TX-1000 Blood Tube Package (GUO)	S4TR refrigerated centrifuge, GUO, 230 V		75015831
	S4TR TX-1000 Blood Tube Package (IVD-MD)	S4TR refrigerated centrifuge, IVD-MD, 230 V		75015827
Floor	X4FR TX-1000 Blood Tube Package (GUO)	X4FR refrigerated floor-standing centrifuge, GUO, 230 V		75015839
	X4FR TX-1000 Blood Tube Package (IVD-MD)	X4FR refrigerated floor-standing centrifuge, IVD-MD, 230 V		75015837
	S4FR TX-1000 Blood Tube Package (GUO)	S4FR refrigerated floor-standing centrifuge, GUO, 230 V		75015823
	S4FR TX-1000 Blood Tube Package (IVD-MD)	S4FR refrigerated floor-standing centrifuge, IVD-MD, 230 V		75015821

Molecular biology applications



Thermo Scientific™ Megafuge™ Centrifuge Packages

Maximum capacities

4x 50 mL	8x 15 mL	24x 2 mL
--------------------	--------------------	--------------------



Image: Megafuge 8 centrifuge



Package	Centrifuge	Accessories included	Cat. No.
Megafuge 8 Ventilated Molecular Biology Package	Megafuge 8 ventilated centrifuge, 230 V	- TX-150 rotor (75005701) – TX-150 round buckets (set of 4) (75005702) – TX-150 ClickSeal biocontainment lids (set of 4) (75005707) – Adapters 50 mL conical (set of 4) (75005736) – Adapters 15 mL conical (set of 4) (75005737) • MicroClick 24 x 2 rotor (75005715)	50147015



Multi-discipline applications



Thermo Scientific™ Megafuge™ Centrifuge Packages

Maximum capacities





Image: Megafuge 8R centrifuge



Package	Centrifuge	Accessories included	Cat. No.
Megafuge 8R Refrigerated Multi-Discipline Package	Megafuge 8R refrigerated centrifuge, 230 V	• TX-150 rotor (75005701) – TX-150 round buckets (set of 4) (75005702) – TX-150 ClickSeal biocontainment lids (set of 4) (75005707) – Adapters 50 mL conical (set of 4) (75005736) – Adapters 15 mL conical (set of 4) (75005737) • HIGHConic III rotor (75005709) • Thermo Scientific™ M10 Microplate Rotor (75005706) – Sealed buckets (set of 2) (7500572) • Thermo Scientific™ MicroClick 24 x 2 Rotor (75005715)	50144443





 Learn more at thermofisher.com/centrifuges

The intended use of the products mentioned in this document varies. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific uses or applications. © 2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **FLY-12488098-EMEA 1125**

thermo scientific