

TW 34XT

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	34XT
LN ₂ Capacity (L)	34
Static Hold Time (Days)	340
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	18.8
Height (in)	26.3
WEIGHT	
Empty (Lbs)	34.75
Full Weight (Lbs)	95.3
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA