

TW 4LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	4LD
LN ₂ Capacity (L)	4
Static Hold Time (Days)	10
Evap. Rate (Liters/Day)	0.40
DIMENSIONS	
Necktube (in/mm)	1.2 / 30.5
Outside Dia. (in/mm)	7.6 / 193
Height (in/mm)	17 / 432
WEIGHT	
Empty (lb/kg)	6.6 / 3.0
Full (lb/kg)	13.7 / 6.2

FEATURES

- 4 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA