

# MTC Bio™

## CULTURE & CENTRIFUGE TUBES



- **Round bottom culture tubes**
- **Conical centrifuge tubes**
- **High clarity polymers**
- **Printed graduations & marking spot**
- **Sterile; RNase & DNase free**



### Culture Tubes

These disposable sterile tubes are made of high quality polypropylene (PP) that will withstand temperatures from -190° to +120°C. They are also available in crystal clear polystyrene (PS) with a more limited range of 0° to 70°C. Caps are made of polyethylene (PE) and have the popular two position design (up/ loose for aerobic work and down/ sealed for anaerobic cultures). Polystyrene tubes may be centrifuged at up to 2000xg. Polypropylene tubes withstand higher forces, up to 4000xg.

### Centrifuge Tubes

Conical 15ml and 50ml centrifuge tubes are made of premium polypropylene for clarity and centrifugal strength. They are supplied sterile, and although they are autoclavable, their polyethylene screw caps are not. Tubes and caps may be frozen at temperatures as low as -45°C. The maximum centrifugal force for both 15ml and 50ml tubes is 9000xg. Caps have a flat top surface for easy marking and labeling. Cat. No. C2605 features a "skirted" bottom that allows the tube to stand upright without a separate test tube rack.



R1030



R1070

**MTC Bio**

Div. of MicroSolv Technology Corp.  
1 Industrial Way West, Bldg. E, Unit D  
Eatontown, NJ 07724 USA  
732-578-1777 • [www.MTCbio.com](http://www.MTCbio.com)

### Ordering Information

#### Culture Tubes (supplied with dual position caps attached)

| Cat. No | Capacity | Dimensions | Material | Packaging Unit    | Units Per Case | Price Per Case | Qty  |
|---------|----------|------------|----------|-------------------|----------------|----------------|------|
| T8301   | 5ml*     | 12x75mm    | PP       | 25/heat seal bag  | 40             | \$159.00       | 1000 |
| T8320   | 5ml*     | 12x75mm    | PP       | 125/heat seal bag | 8              | \$99.00        | 1000 |
| T8300   | 5ml*     | 12x75mm    | PS       | 125/heat seal bag | 8              | \$99.00        | 1000 |
| T8235   | 14ml*    | 17x100mm   | PP       | 25/notched bag    | 40             | \$209.00       | 1000 |
| T8260   | 14ml*    | 17x100mm   | PP       | 1/sterile peel    | 500            | \$129.00       | 500  |
| T8225   | 14ml*    | 17x100mm   | PS       | 25/notched bag    | 40             | \$215.00       | 1000 |
| T8250   | 14ml*    | 17x100mm   | PS       | 1/sterile peel    | 500            | \$129.00       | 500  |

\* 12x75 tubes are graduated to 4ml. 17x100 tubes are graduated to 12ml

#### Centrifuge Tubes (supplied with screw caps attached)

| Cat. No | Capacity                                                                            | Dimensions | Material | Packaging Unit    | Units Per Case | Price Per Case/500 |
|---------|-------------------------------------------------------------------------------------|------------|----------|-------------------|----------------|--------------------|
| C2600   | 15ml                                                                                | 17x118mm   | PP       | 50 in foam rack   | 10             | \$129.00           |
| C2601   | 15ml                                                                                | 17x118mm   | PP       | 50 in notched bag | 10             | \$119.00           |
| C2602   | 50ml                                                                                | 29x115mm   | PP       | 25 in foam rack   | 20             | \$149.00           |
| C2603   | 50ml                                                                                | 29x115mm   | PP       | 25 in notched bag | 20             | \$139.00           |
| C2605*  | 50ml                                                                                | 29x115mm   | PP       | 25 in notched bag | 20             | \$159.00           |
| R1070   | Rack for 30x15ml tubes, pack of 5                                                   |            |          |                   |                | \$50.00            |
| R1030   | 4-way rack for 4 x 50ml, 12 x 15ml, 32 x 1.5ml, 32 x 0.5ml tubes, rainbow pack of 5 |            |          |                   |                | \$50.00            |

\* No. 2605 is a self-standing, sterile tube